



A Step by Step of Ultrasonic Microscope Design

By:
Agus Indra Gunawan

WEBINAR SERIES
(seri-8, 02 September 2020)

1

1

Bidang Keilmuan

- **Ultrasonic and Non Destructive Measurement**
Strength of Clay, Concrete, Biodiesel, Oil, Salinity
- **Transesterification (Converting Jelantah to Biodiesel)**
Transesterification machine, using ultrasonic as a vibration generator
- **Aquaculture Engineering Application**
Water quality, shrimp seed counter, DO sensor, salinity sensor, feeding management system, etc.


<https://aigunawan.blogspot.com/>






2

What is Ultrasonic ??

- Gelombang mekanik
- Getaran dari partikel











3

Ultrasonic today







- In the nature → Bat, Dolphin
- Robot obstacle avoidance
- Ultrasonography
- Car parking guidance
- Machine Cleaner
- Level Detector
- Crack material



[1]. <https://www.honda-el.co.jp/en/medical/HS-2600.html>



4



5

Perbedaan Optical dan Acoustical Observation

Need to stain to obtain
a good contrast

Kidney specimen 800 x 800 pixel
(H&E stained)

After staining, the contrast of specimen
is much better than unstained one.

- High resolution.
- It is not easy to obtain a good contrast.
- Staining process.
- It needs several hours to days to complete the staining.
- **Severe staining could damage the biological function of the specimen.**

6




Perbedaan Optical dan Acoustical Observation

- Rapid observation can be realized by skipping the staining process.
- The contrast represents visco-elastic property.
- In vitro and in vivo (living specimen) observation.
- Contaminant free (some cases)
- High frequency is required for high resolution.



7

Persiapan Sampel




Speed of sound is 2340 m/s
Ac. Impd is 2.46 MNs/m³




8

Acoustical and Optical Result




9

Normal gastric structure




10

Normal gastric mucosa and submucosa




11


Glia and Glioma Observation




12




13




14




15




16




17




18




19




20




21



22



23




24

Sound Field Analysis

Improvement:

- Upgrading spatial resolution by replacing 16° transducer with 22° transducer
- All aspects including oblique incidence, pressure wave and shear wave are also considered
- Calculate sound field by utilizing Fourier analysis in k -space




25

Sound Field Analysis

$$\phi(x, y, z) = \phi_0 \frac{e^{-j \cdot (k \cdot l)}}{|l|} da$$

$$\mathbf{k} = k_x + k_y + k_z$$

$$k = \frac{c}{\omega}$$




26

Sound Field Analysis

$$\phi(x,y,z) = \phi_0 \frac{e^{-j \cdot (k \cdot l)}}{|u|} da$$

$$\mathbf{k} = k_x + k_y + k_z$$

$$k = \frac{c}{\omega}$$





27

27


Sound Field Analysis

Sound field distribution (2D)




28

28




29



30

Oblique Incidence Wave



31


Oblique Incidence Wave (TF)

$$\begin{bmatrix} k_{s1} & k_{p1} & k_{s2} & -k_{p2} \\ -k_{s1} & k_{p1} & k_{s2} & k_{p2} \\ \rho_1 & -\rho_1 & \rho_2 & \rho_2 \\ -\rho_1 & -\frac{\mu_1 \rho_1}{\lambda_1 + 2\mu_1} & -\rho_2 & -\frac{\mu_2 \rho_2}{\lambda_2 + 2\mu_2} \end{bmatrix} \begin{bmatrix} \cos \theta_{sr} & \sin \theta_{pr} & \cos \theta_{st} & \sin \theta_{pt} \\ \sin \theta_{sr} & \cos \theta_{pr} & \sin \theta_{st} & \cos \theta_{pt} \\ \sin 2\theta_{sr} & \cos 2\theta_{sr} & \sin 2\theta_{st} & \cos 2\theta_{st} \\ \cos 2\theta_{sr} & \sin 2\theta_{pr} & \cos 2\theta_{st} & \sin 2\theta_{pt} \end{bmatrix} \begin{bmatrix} A_{sr} \\ A_{pr} \\ A_{st} \\ A_{pt} \end{bmatrix} =$$


$$\begin{bmatrix} k_{s1} & -k_{p1} \\ k_{s1} & k_{p1} \\ \rho_1 & \rho_1 \\ \rho_1 & -\frac{\mu_1 \rho_1}{\lambda_1 + 2\mu_1} \end{bmatrix} \begin{bmatrix} \cos \theta_{sr} & \sin \theta_{pr} \\ \sin \theta_{sr} & \cos \theta_{pr} \\ \sin 2\theta_{sr} & \cos 2\theta_{sr} \\ \cos 2\theta_{sr} & \sin 2\theta_{pr} \end{bmatrix} \begin{bmatrix} A_{si} \\ A_{pi} \end{bmatrix}$$



32



33



34

Reflection Intensity

TRANSDUCER

```

graph TD
    Start([start]) --> Determine[Determine 2-D k-space (k_x, k_y)]
    Determine --> Create[Create Transfer Function]
    Determine --> FFT[2-Dimensional Fourier Transform (FFT_{rela})]
    Create --> Transfer[2-D Transfer Function]
    FFT --> Shifted["Shifted Fourier Transform both two medium  
FFT_0 = e^{(-j*k_z*d)} * FFT_{field}"]
    Shifted --> Multiply["Multiply & sum operation"]
    Transfer --> Multiply
    Multiply --> Inverse["2-Dimensional Invers Fourier Transform (IFFT_{rela})"]
    Inverse --> Plot[Plot the graph]
    Plot --> Stop([stop])
    
```

35

35

Reflection Intensity


```

graph TD
    Start([start]) --> Determine[Determine 2-D k-space (k_x, k_y)]
    Determine --> Create[Create Transfer Function]
    Determine --> FFT[2-Dimensional Fourier Transform (FFT_{rela})]
    Create --> Transfer[2-D Transfer Function]
    FFT --> Shifted["Shifted Fourier Transform both two medium  
FFT_0 = e^{(-j*k_z*d)} * FFT_{field}"]
    Shifted --> Multiply["Multiply & sum operation"]
    Transfer --> Multiply
    Multiply --> Inverse["2-Dimensional Invers Fourier Transform (IFFT_{rela})"]
    Inverse --> Plot[Plot the graph]
    Plot --> Stop([stop])
    
```

36

36

Acoustic Impedance Profile of a Glia in Enlarged Scale



This figure was 250 x 250 μm, covered by 200 x 200 pixels.

Cell Motility, defined by Mitchison and Cramer [103], Raftopoulou and Hall [105].


Classified into 3 areas:

1. Protrusion area → microtubules
2. Cell body → the thickest region containing nucleus (rich of actin filament and microtubules).
3. Intermediate area → lay in between cell body and protrusion area (rich of actin filament and microtubules).



37

A Trace in Acoustic Impedance Along a Line



Either Buffer Liquid or Pure water can be used as a reference



38

Collaborative research and Cooperation are
Welcome

**Terima kasih
atas perhatiannya**

